

## Urine Collection Options for FedEx Pick up:

Collection Media	Test	Viable for testing from time of collection
Red/Yellow Marble Top BD Vacutainer w/ Preservative	Urinalysis Only	72 Hours
Grey top with preservative (C&S Transfer Straw Kit)	Urine Culture	48 Hours
Urine Collection Kit	Urinalysis and Urine Culture	See 2 items above (marble and grey top tubes)
Yellow Screw top Boritex Tube	Urinalysis and Urine Culture	72 Hours



**Red & Yellow Marble BD Vacutainer/UA Preservative tube – Urinalysis Only.**

The sample is viable for testing 72 hours after time of collection.



**C&S Transfer Straw Kit – Urine Culture Only.** The sample is viable for testing 48 hours after time of collection.



**Urine Collection Kit – Urinalysis and Urine Culture collection.**

The red/yellow urinalysis tube included in the urine collection kit is viable for testing 72 hours from collection. The grey top tube is for urine culture and is viable for 48 hours from the time of collection.



**Boritex Tube – Urinalysis and Urine Culture.**

The sample is viable for testing 72 hours after time of collection.